SELF-ETCH DENTIN ENAMEL BONDING LIGHT CURE

INSTRUCTIONS FOR USE

English
Español
Français

 CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd.-1st.
P.C. 08025, Barcelona, Spain.

Rev. Rev."0" / 21 Jul 25

SELF-ETCH DENTIN ENAMEL BONDING LIGHT CURE

ENGLISH

PRODUCT DESCRIPTION: MEDENTAL Bonding Agents are monomer-based composite materials including Bis-GMA, UDMA, HEMA and TEGDMA that are used for bonding dental materials such as composite resins, ceramics or cements to enamel or dentin.

CLINICAL BENEFITS:

- Improved aesthetics.
- Conservation of the restoration.
- They reduce sensitivity.
- Simple repairs.
- Durability.
- It does not give off harmful components; it is biocompatible with the body.

PERFORMANCE CHARACTERISTICS:

- Good adhesion.
- Marginal integrity.
- Holding force.
- Compatibility with other restorative materials.
- For better results it is recommended to use halogen or LED lamp with a wavelength between 400-500nm that activates the photoinitiators of the bonding. For better results is recommended to use halogen or LED lamp with a wavelength between 400-500nm that activates the photoinitiators of the bonding.

INTENDED USE: To enhance the bonding strength between the tooth enamel and the composite, and to seal the marginal edges.

INDICATIONS:

Bonding agents are indicated for all kind of restorations with composites and compomers, cements, ceramics and endodontic treatments.

INTENDED USER: Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must make restorations with Medental Self-Etch Dentin Enamel Bonding Light Cure.

PRECAUTIONS:

Some patients can present allergies to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away from the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

No requires refrigeration.

Only for professional dental use.

Room temperature before use.

Not require cold chains to transport.

Do not use it after the date of expiry.

CONTRAINDICATIONS:

-The material described must not be used when the application techniques are not possible.

- Bonding Agents are subject to irritation in soft tissues and cells in direct contact when the monomers included in their formulation are in direct contact with these tissues, causing irritation and even pulp death if they are not 100% photocured.

- It is contraindicated to place bonding agent in any presentation on materials containing eugenol since these inhibit polymerization.

- Avoid contamination with saliva or blood.

- Monitor the moisture content of the dentin and its compatibility with the material to be used.

COMPOSITION: PMGDM in acetone; HEMA; Bidistilled water; Diurethane; Photo initiator, Stabilizer.

INSTRUCTIONS FOR USE:

A. Tooth preparation: Proper isolation of the tooth and the area to be worked on must be carried out using clamps and a rubber dam.

B. Enamel etching (For bonding materials requiring etching): Apply 35% phosphoric acid to the enamel and

dentin for 20 to 30 seconds on enamel and 10 to 15 seconds on dentin. Rinse with water and ensure the enamel has a grayish appearance, keeping the tissue moist but not drying out. If necessary, repeat the process. THIS STEP SHOULD BE OMITTED WHEN USING SELF-ETCHING BONDING MATERIAL. The application of etching acid is done for the purpose of opening dentin tubules and creating retention.

C. Place a drop of dentin enamel bonding agent on a glass slab and apply a thin layer to cover the dentin and etched enamel.

D. Apply a thin layer of enamel adhesive with a brush and use a gentle stream of moisture- or oil-free air to evenly disperse the thin film of adhesive over the restoration area. Then light-cure for 10–20 seconds. This will create a bond between the tooth and the restoration.

E. With a brush or cotton pellet take liquid from the slab and apply a fine layer of the over the dentin and etched enamel.

F. When the application of the bonding is finished, cure it with light cure lamp (LED 800mw/cm²) for 20 seconds.
G. Do not touch the remaining wet film of the bonding layer, and then proceed to make the restoration with the selected material.

CONSERVATION:

- Close the bottle immediately after to dispense the product.
- To prevent deterioration of material, keep bottle tightly closed when not in use.
- Do not expose materials to elevated temperature or intense light.
- The system is designed to be used at room temperature above 24°C or 75°F.

WARNINGS and PRECAUTIONS:

- Read carefully and follow all instructions
- Before the use of these products always use rubber dam gloves and protective glasses.
- The etching solution contains a 37% phosphoric acid solution, which can be harmful to skin and eyes.

-In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor.

-In case of contact with skin wash the area with soap and water.

-The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance.

-Ask for patient to rinse his mouth immediately after treatment.

-The product is designed to be used at room temperature.

-It is recommended to cover the polymerization time indicated by the manufacturer, according to the layer thickness, because a non-polymerized bonding or resin can cause irritation or allergic reactions.

- On case of Self-etch is not necessary to apply acid etchant.

- Not all Universal bondings are indicated to all types of restorative materials.

- The Etch and rinse adhesives require careful technique and may be sensitive to the level of dentin moisture.

- The bond strength of Self-Etch bonding may be slightly lower compared to Etch and Rinse systems.

- Dentin thickness positively affected biocompatibility.

WARRANTY:

If the use of this product presents some anomaly such as:

- It does not fulfill the intended purpose,
- It doesn't return the functionality of the tooth,
- It doesn't have the properties it refers to,
- It has damage to packaging or manufacturing,
- The instructions for use are misunderstood,
- Any secondary reaction or anomaly:

Please send detailed report to: Conamco, S.A. de C.V. email: atencionclientes@conamcosacv.com.mx

Warranty: MEDENTAL and/or CONAMCO S.A. de C.V. warranty the product against any fault of manufacture of the material. MEDENTAL and/or CONAMCO S.A. de C.V. do not take responsibility for the use of the product

outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective products.

The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at:

<https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS:

1. Use New and clean nitrile gloves.
2. The Bonding and resin must be cured in a safe container for resin, either by exposing it to the sun or with a light cure lamp.
3. Once Cured, it can be discarded with solid residues.

NOTE: Liquid or partially cured resin or bonding must be disposed of in accordance with the regulations.

FEATURES/CLAIMS:

- Similar biocompatibility.
- Higher antimicrobial activity against streptococcus mutans.
- Good shear bond straight values.
- More resistance and less operator sensitive.
- Good adhesion.
- Prior use of adhesives to occlude dentinal tubules could reduce the discoloration induced by MTA.
- Effective reducing dentin hypersensitivity.

To take good care of the restoration, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and use dental protector if the patient has Bruxism.

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- Veuillez envoyer un rapport détaillé à : Conamco, S.A. de C.V., email :