GLASS IONOMER CEMENT TYPE II

RESTORATIVE ART (Atraumatic Restoration Technique)

INSTRUCTIONS FOR USE

English
Español
Français

 CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd. -1st.
P.C. 08025, Barcelona, Spain.

Rev: "0" / 21 JUL 2025

GLASS IONOMER CEMENT TYPE II (RESTORATIVE ART)

ENGLISH

PRODUCT DESCRIPTION: Medental Glass Ionomer cements (also called polyalkenoate cements) are a category of the cements group used by Dentists. They are made by diluted polyacrylic acid and fine glass powder without lead. Polyacrylics form with the glass a strong matrix chemically compatible with dental composite, metal, and dentine. They do not shrink, which is why these cements are successfully used as luting, liners or restorative depending on the thickness of the polyacrylic acid liquid in the mixture with the glass. Glass Ionomer is a "long term cement" of use from one week to four months but also, they can be used as definitive cements.

CLINICAL BENEFITS:

- Minimally invasive technique
- Caries removal with hand instruments
- Restoration with high viscosity GIC

PERFORMANCE CHARACTERISTICS:

- Chemical adhesion to enamel and dentin.
- Menor tasa de descalcificación
- Chemical bonding without etching.
- Fluoride release.

COLOR: Universal Tooth shade.

INTENDED USE:

For the atraumatic restorative techniques. Permanent restoration of deciduous teeth, class III & V, and cervical erosion. Long term temporary restorative for Class I and II.

INDICATIONS:

- Children
- Patients with limited access to dental care
- Environments without electricity or rotating equipment

INTENDED USER: Dental professionals only with an academic degree, experience, and training for use

dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with primary or permanent dentition, who, for different reasons, must be make restorations with Medental Glass Ionomer Type II Restorative ART.

PRECAUTIONS: Some patients may present allergy to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

Does not require refrigeration.

Only for professional dental use.

Room temperature before use. It does not require cold chain to transport.

Do not use after date of expiry.

CONTRAINDICATIONS:

-The material must not be used when the application techniques described are not possible.

- Very deep lesions with pulp involvement

- Large cavities in loading surfaces

- Presence of abscess, fistula, or spontaneous pain

- Teeth with mobility or irreversible pulp pathology

-If there is direct exposure of the Pulp is no recommended to use the glass ionomer alone.

-Do not use if efficient isolation cannot be achieved, due to its sensitivity to humidity.

COMPOSITION: Alumino-silicated glass, polyacrylic acid, varnish.

INSTRUCTIONS FOR USE:

A. POWDER AND LIQUID DISPENSING:

1. The standard powder to liquid ratio is 1.2g of powder and 0.4ml of liquid.

2. For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
3. Hold the liquid bottle vertically and squeeze gently.
4. Close bottles immediately after use.
- B. **MIXING:** Dispense powder and liquid onto the pad. Using the plastic spatula. Divide the powder into 2 equal parts. Mix the first portion with all the liquid for 10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-20 seconds (total time within 30 seconds).
- C. **PRODUCT PLACEMENT:** Place the mixture in the prepared cavity (the cavity has to be dry and clean). When all the restoration is met, let the product lose its brightness, then put your finger over the restoration pushing inside the cavity between one or two minutes. When 6 minutes have elapsed, cover the restoration with layer of varnish (The varnish avoids water absorption of the product. The product needs at least 24 hours to meet the maximum strength)

CONSERVATION:

- Close the package immediately after to dispense the product.
- To prevent deterioration of material, keep bottle and jar tightly closed when not in use.
- Replace the top of the bottle immediately after use to prevent contamination.
- Do not store in places exposed to humidity.

WARNINGS and PRECAUTIONS: In case of any change to the original function of the product, stop using it and contact the manufacturer.

-The Glass Ionomer Type II Restorative ART is not irritable to the skin.

- Read carefully and follow all instructions
 - Before the use of these products always use rubber dam gloves and protective glasses.
 - In case of contact with the eye, wash immediately with abundant water for at least 15 minutes and if the irritation persists consult a doctor.
 - In case of contact with skin wash the area with soap and water.
 - The product should not be ingested or aspirated. If accidentally swallowing induces vomit immediately and get medical attention. In case of inhalation, clean the nose with water and seek medical assistance.
 - Ask for patient to rinse his mouth immediately after treatment.
 - The product is designed to be used at room temperature.
 - The remains of cement in the dentin can interfere with adhesion if a resin is placed on the same place.
 - No cement meets all clinical needs, it is the responsibility of the professional the choice of cement that meets the clinical requirements as the case may be.
- WARRANTY:** If the use of this product on a patient or anyone else causes some anomaly, please send detailed report to: Conamco, SA de CV email:

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Warranty: MEDENTAL INTERNATIONAL and/or CONAMCO SA de CV warranty the product against any fault of manufacturing of the material. MEDENTAL INTERNATIONAL AND/OR CONAMCO SA de CV does not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product. The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at:

<https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS: New, clean nitrile gloves should be worn.

- 1) *Removal of excess.*

- Fresh (before complete setting): Remove with a probe or spatula. recommends cleaning immediately.
- Setting: Using cement chisels, cures or metal excavators with controlled movements, pass waxed dental floss in interproximal areas, pulling downwards to avoid fracturing restorations.

2) *Removal under restorations or crowns.*

- It is not possible to remove it without removing the restoration; in these cases, the crown must be sectioned with a carbide/diamond bur and the tooth cleaned manually or with ultrasound.

*In both cases, it is recommended to protect soft tissues, work with good lighting, and use a magnifying glass to identify cement residues.

- 3) Dispose of with solid waste in accordance with the regulations of each country.

FEATURES CLAIMS: Good mechanical resistance.

- Minimally invasive.
 - Anticariogenic.
 - Fluoride release.
 - Good chemical adhesion to enamel and dentin.
 - No etching required.
 - Thermal and chemical isolation.
 - Biocompatible.
 - Low toxicity.
 - Prevents descaling
 - High success rate in Class I cavities in primary dentition
 - Moderate in class II jar tightly closed when not in use.
 - Excellent in public health programs.
- To take good care of the restoration, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and use dental protector if the patient has Bruxism.

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