GLASS IONOMER CEMENT TYPE I

LUTING

INSTRUCTIONS FOR USE

English
Español
Français

 CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd.-1st.
P.C. 08025, Barcelona, Spain.

Rev: "0" / 21 JUL 2025

GLASS IONOMER CEMENT TYPE I (LUTING)

ENGLISH

PRODUCT DESCRIPTION: Medental Glass Ionomer cements (also called polyalkenoate cements) are a category of the cements group used by Dentists. They are made by diluted polyacrylic acid and fine glass powder without lead. Polyacrylics form with the glass a strong matrix chemically compatible with dental composite, metal, and dentine. They do not shrink, which is why these cements are successfully used as luting, liners or restorative depending on the thickness of the polyacrylic acid liquid in the mixture with the glass. Glass Ionomer is a "long term cement" of use from one week to four months but also, they can be used as definitive cements.

CLINICAL BENEFITS:

- Symptoms relief (pain and discomfort)
- Esthetic (Universal tooth shade)
- Radiolucent
- Rebuild and return chewing functionality
- It does not give off harmful components; it is biocompatible with the body.

PERFORMANCE CHARACTERISTICS:

- High adhesion
- Elastic Model
- Resistance
- Long lasting bonding and strength integrity

COLOR: Universal Tooth shade.

INTENDED USE:

Restorations in class V and low charge class III or small class I cervical lesions, temporary or intermediate restorations as a preventive (sealant) in patients with high risk of caries.

INDICATIONS:

Cementation of cavities Class V and III; inlays, onlays in children and adults.

INTENDED USER: Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with primary or permanent dentition, who, for different reasons, must be make restorations with Medental Glass Ionomer Type I Luting Cement.

PRECAUTIONS: Some patients may present allergy to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

Does not require refrigeration.

Only for professional dental use.

Room temperature before use. It does not require cold chain to transport.

Do not use after date of expiry.

CONTRAINDICATIONS:

-The material must not be used when the application techniques described are not possible.

-Do not use in extensive class II high occlusal load cavities.

-Not recommended for demanding aesthetic restorations.

-If there is direct exposure of the Pulp is no recommended to use the glass ionomer alone.

-Do not use if efficient isolation cannot be achieved, due to its sensitivity to humidity.

COMPOSITION: Alumino-silicated glass, polyacrylic acid.

INSTRUCTIONS FOR USE:

- a. The standard powder to liquid ratio is 0.85g of powder and 0.45g of liquid.

- b. For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
- c. Hold the liquid bottle vertically and squeeze gently.
- d. Close bottles immediately after use.
- e. **MIXING:** Dispense powder and liquid onto the pad. Using the plastic spatula. Divide the powder into 2 equal parts. Mix the first portion with all the liquid for 10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-20 seconds (total time within 30 seconds).

CONSERVATION:

-Close the package immediately after to dispense the product.

-To prevent deterioration of material, keep bottle and jar tightly closed when not in use.

-Replace the top of the bottle immediately after use to prevent contamination.

-Do not store in places exposed to humidity.

WARNINGS and PRECAUTIONS: In case of any change to the original function of the product, stop using it and contact the manufacturer.

-The Glass Ionomer Type I Luting Cement is not irritable to the skin.

-Read carefully and follow all instructions

-Before the use of these products always use rubber dam gloves and protective glasses.

-In case of contact with the eye, wash immediately with abundant water for at least 15 minutes and if the irritation persists consult a doctor.

-In case of contact with skin wash the area with soap and water.

-The product should not be ingested or aspirated. If accidentally swallowing induces vomit immediately and get medical attention. In case of inhalation, clean the nose with water and seek medical assistance.

-Ask for patient to rinse his mouth immediately after treatment.

-The product is designed to be used at room temperature.

- The remains of cement in the dentin can interfere with adhesion if a resin is placed on the same place.

- No cement meets all clinical needs, it is the responsibility of the professional the choice of cement that meets the clinical requirements as the case may be.

WARRANTY: If the use of this product on a patient or anyone else causes some anomaly, please send detailed report to: Conamco, SA de CV email:

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Warranty: MEDENTAL INTERNATIONAL and/or CONAMCO SA de CV warranty the product against any fault of manufacturing of the material. MEDENTAL INTERNATIONAL AND/OR CONAMCO SA de CV does not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product. The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at: <https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS: New, clean nitrile gloves should be worn.

1) *Removal of excess.*

- Fresh (before complete setting): Remove with a probe or spatula. recommends cleaning immediately.

- Setting: Using cement chisels, cures or metal excavators with controlled movements, pass waxed dental floss in interproximal areas, pulling downwards to avoid fracturing restorations.

2) *Removal under restorations or crowns.*

- It is not possible to remove it without removing the restoration; in these cases, the crown must be sectioned with a carbide/diamond bur and the tooth cleaned manually or with ultrasound.

*In both cases, it is recommended to protect soft tissues, work with good lighting, and use a magnifying glass to identify cement residues.

3) Dispose of with solid waste in accordance with the regulations of each country.

FEATURES CLAIMS: Good mechanical resistance.

-Minimally invasive.

-Anticariogenic.

-Fluoride release.

-Good chemical adhesion to enamel and dentin.

-No etching required.

-Useful in deep cavities without pulp exposure.

-Thermal and chemical isolation.

-Biocompatible.

-Low toxicity.

To take good care of the restoration, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and use dental protector if the patient has Bruxism.

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