GLASS IONOMER CEMENT TYPE I ORTHODONTIC

INSTRUCTIONS FOR USE

English
Español
Français

 CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd.-1st.
P.C. 08025, Barcelona, Spain.

Rev: "0" / 21 JUL 2025

GLASS IONOMER CEMENT TYPE I (ORTHODONTIC)

ENGLISH

PRODUCT DESCRIPTION: Medental Glass Ionomer cements (also called polyalkenoate cements) are a category of the cements group used by Dentists. They are made by diluted polyacrylic acid and fine glass powder without lead. Polyacrylics form with the glass a strong matrix chemically compatible with dental composite, metal, and dentine. They do not shrink, which is why these cements are successfully used as luting, liners or restorative depending on the thickness of the polyacrylic acid liquid in the mixture with the glass. Glass Ionomer is a "long term cement" of use from one week to four months but also, they can be used as definitive cements.

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- Cementation of orthodontic bands

INDICATIONS:

Cementation of orthodontic bands with caries prevention.

INTENDED USER: Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with primary or permanent dentition, who, for different reasons, must be make restorations with Medental Glass Ionomer Type I Orthodontic.

PRECAUTIONS: Some patients may present allergy to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

Does not require refrigeration.

Only for professional dental use.

Room temperature before use. It does not require cold chain to transport.

Do not use after date of expiry.

CONTRAINDICATIONS:

-The material must not be used when the application techniques described are not possible.

- Situations where high immediate adhesive strength is required (preferably resin)
- Surfaces with poor band adaptation
- Patients with very high occlusal forces
- Slower setting or lower adhesion compared to resins

COMPOSITION: Alumino-silicated glass, polyacrylic acid.

INSTRUCTIONS FOR USE:

- Prior to mixing, shake powder bottle several times to fluff.
- Dispense 1 level scoop of Orthodontic powder, or multiples there of on to mixing pad.
- Dispense 2 droops of liquid onto mixing pad for each 1 scoop of powder.
- Using a suitable spatula, incorporate the Orthodontic powder into the liquid and thoroughly mix with a flattering action for 20 seconds.

e. The Orthodontic Cement has a working time of 3 minutes.

f. **APPLICATION:**

1. Clean teeth to be banded with pumice or non-oil containing prophylactic paste and rinse thoroughly.
2. Dry thoroughly and provide suitable moisture control to isolate area and to prevent moisture contamination during the initial setting.
3. Apply Orthodontic Cement onto all tooth bearing surfaces of the band.
4. Immediately seat band onto tooth.
5. The Orthodontic cement can be easily wiped away with a cotton roll during the gel stage or removed with a scaler after it has set.

CONSERVATION:

-Close the package immediately after to dispense the product.

-To prevent deterioration of material, keep bottle and jar tightly closed when not in use.

-Replace the top of the bottle immediately after use to prevent contamination.

-Do not store in places exposed to humidity.

WARNINGS and PRECAUTIONS: In case of any change to the original function of the product, stop using it and contact the manufacturer.

-The Glass Ionomer Type I Luting Cement is not irritable to the skin.

-Read carefully and follow all instructions

-Before the use of these products always use rubber dam gloves and protective glasses.

-In case of contact with the eye, wash immediately with abundant water for at least 15 minutes and if the irritation persists consult a doctor.

-In case of contact with skin wash the area with soap and water.

-The product should not be ingested or aspirated. If accidentally swallowing induces vomit immediately and

get medical attention. In case of inhalation, clean the nose with water and seek medical assistance.

-Ask for patient to rinse his mouth immediately after treatment.

-The product is designed to be used at room temperature.

-The remains of cement in the dentin can interfere with adhesion if a resin is placed on the same place.

-No cement meets all clinical needs, it is the responsibility of the professional the choice of cement that meets the clinical requirements as the case may be.

WARRANTY: If the use of this product on a patient or anyone else causes some anomaly, please send detailed report to: Conamco, SA de CV email:

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Warranty: MEDENTAL INTERNATIONAL and/or CONAMCO SA de CV warranty the product against any fault of manufacturing of the material. MEDENTAL INTERNATIONAL AND/OR CONAMCO SA de CV does not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product. The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at:

<https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS: New, clean nitrile gloves should be worn.

1) *Removal of excess.*

- Fresh (before complete setting): Remove with a probe or spatula. recommends cleaning immediately.

- Setting: Using cement chisels, cures or metal excavators with controlled movements, pass waxed dental floss in interproximal areas, pulling downwards to avoid fracturing restorations.

2) *Removal under restorations or crowns.*

- It is not possible to remove it without removing the restoration; in these cases, the crown must be sectioned

with a carbide/diamond bur and the tooth cleaned manually or with ultrasound.

*In both cases, it is recommended to protect soft tissues, work with good lighting, and use a magnifying glass to identify cement residues.

3) Dispose of with solid waste in accordance with the regulations of each country.

FEATURES CLAIMS: Good mechanical resistance.

-Minimally invasive.

-Anticariogenic.

-Fluoride release.

-Good chemical adhesion to enamel.

-No etching required.

-Thermal and chemical isolation.

-Biocompatible.

-Low toxicity.

-Prevents descaling

To take good care of the restoration, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and use dental protector if the patient has Bruxism.

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