PROVISIONAL CEMENT

EUGENOL FREE

INSTRUCTIONS FOR USE

English
Español
Français

 CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd.-1st.
P.C. 08025, Barcelona, Spain.

Rev: "0" / 21 JUL 2025

PROVISIONAL CEMENT EUGENOL FREE

ENGLISH

PRODUCT DESCRIPTION: MEDENTAL Temporary Cements Eugenol Free is recommended for use in a short time and when the definitive restoration is a light cure cement or resin. Most prosthesis can be removed simply by twisting or applying a small force with a crown remover. The provisional cement is a temporary filling for self-curing under moisture.

CLINICAL BENEFITS:

- Avoid allergic reactions or sensitivity.
- Easy removal.
- Compatibility with modern restorations.
- Good temporary sealing.
- Rapid initial hardening (60-80 sec.)
- It does not stain temporary acrylic teeth.

PERFORMANCE CHARACTERISTICS:

• The permanent cementation doesn't have to exceed the estimated time to use non-eugenol cement. • They have lower retention than other cements or resin cements, especially when Retrievability is desired. • Mechanical and chemo-mechanical cleaning protocols of non-eugenol cements can influence restoration bond strength. Clinically relevant sealing for provisional restorations in short-term simulations. • Microleakage in Zinc Oxide Non-eugenol is more in comparison with others cements. Adequate mechanical resistance. More homogeneous texture. Reduced solubility. Moderate short-term compressive strength.

COLOR: Dentin color (light beige/cream).

INTENDED USE:

For occlusal temporary seals. For temporary filling after endodontic treatments.

INDICATIONS:

- As a temporary cement in cases of adhesive indirect restorations such as inlays, veneers and crowns.

- When it is necessary to fully preserve the adhesive potential of dentin and the internal restoration substrate.

- In cementation of intraradicular posts of fiberglass or carbon.

- On cast metal posts.

- In orthodontic cases where temporary crowns are required for several months.

- In patients with sensitivity or allergic reactions to eugenol.

INTENDED USER: Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must be make restorations with Medental Provisional Cement (Eugenol Free).

PRECAUTIONS: Some patients may present allergy to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

Does not require refrigeration.

Only for professional dental use.

Room temperature before use. It does not require cold chain to transport.

Do not use after date of expiry.

CONTRAINDICATIONS:

• Do not use if the tooth is very sensitive post-operatively. • If a sedative effect is sought in a preparation that causes discomfort. • If it is required that the cement be removed more easily. • If long-term permanent cementation is required. • If pulp irritation exists in deep cavities. • Postoperative discomfort in

deep areas. • Avoid in crowns, bridges, posts and definitive or permanent restorations. • Avoid use in highly sensitive pediatric patients. • If there is an allergy to any of the components of the formula.

COMPOSITION: Calcium Sulfate, Zinc oxide, Zinc stearate, cotton oil.

INSTRUCTIONS FOR USE:

A) Using an instrument, apply necessary amount in the wet cavity and cover the bottle immediately.

B) The material hardens after a few minutes.

C) Avoid occlusion loads until 2 hours of application.

D) Deep cavities do not need to be filled to the bottom. When sealing a cure, be careful not to exert pressure when applying.

E) To remove the temporary obturation, use rotary instruments.

CONSERVATION:

- To prevent deterioration of material, keep the jar tightly closed when not in use.

-Due to its sensitivity to humidity, cover the bottle immediately after to dispense the product, otherwise it would harden due to the entry of moisture and probably the material could not be extracted again.

WARNINGS and PRECAUTIONS: In case of any change on the original function of the product, stop using it and contact the manufacturer.

- Read carefully and follow all instructions.

- It may cause mild pulp irritation or inflammation immediately after its application.

- It can cause local allergies and contact dermatitis. Its clinical use is relatively safe if handled properly.

- It is recommended to thoroughly clean the remnants before using the final cement or resin.

- Should not be used as a permanent cement.

- Should not be used if a strong bond is required.

- In deep cavities, place a protective intermediate base to avoid irritation.

- Before the use of these products always use rubber dam gloves and protective glasses.

- In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor.

- In case of contact with skin wash the area with soap and water.

- The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance.

- Ask for patient to rinse his mouth immediately after treatment.

- The product is designed to be used at room temperature.

WARRANTY: If the use of this product on a patient or anyone else causes some anomaly, please send detailed report to: Conamco, SA de CV email:

atencionaclientes@conamcosacv.com.mx

Warranty: MEDENTAL INTERNATIONAL and/or CONAMCO SA de CV warranty the product against any fault of manufacturing of the material. MEDENTAL INTERNATIONAL AND/OR CONAMCO SA de CV does not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product. The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at: <https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS: New, clean nitrile gloves should be worn.

1) *Removal of excess.*

- Fresh (before complete setting): Remove with a probe or spatula. recommends cleaning immediately.

- Setting: Using cement chisels, curettes or metal excavators with controlled movements, pass waxed dental floss in interproximal areas, pulling downwards to avoid fracturing restorations.

2) *Removal under restorations or crowns.*

- It is not possible to remove it without removing the restoration; in these cases, the crown must be sectioned with a carbide/diamond bur and the tooth cleaned manually or with ultrasound.

*In both cases, it is recommended to protect soft tissues, work with good lighting, and use a magnifying glass to identify cement residues.

3) Dispose of with solid waste in accordance with the regulations of each country.

FEATURES CLAIMS: Zinc Oxide Free-Eugenol based cements are used as temporary cements short-term. Eugenol-free temporary cements are effective in the clinical setting in the short term. Since they are not used as permanent solutions, they are easy to remove, but their performance is limited to short-term, so it is recommended not to leave too much time between temporary and permanent restoration. It is necessary to ensure that no remaining funds to ensure the adhesion of the final restorative material. To take good care of the restoration, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and use dental protector if the patient has Bruxism.

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