TEMPORARY CEMENT

EUGENOL FREE

INSTRUCTIONS FOR USE

English
Español
Français

 CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd.-1st.
P.C. 08025, Barcelona, Spain.

Rev: "0" / 21 JUL 2025

TEMPORARY CEMENT EUGENOL FREE

ENGLISH

PRODUCT DESCRIPTION: MEDENTAL Temporary Cements Eugenol Free is recommended for use in a short time and when the definitive restoration is a light cure cement or resin. Most prosthesis can be removed simply by twisting or applying a small force with a crown remover. The provisional cement is a temporary filling for self-curing under moisture.

The Medental temporary cement eugenol free is a base of Zinc oxide, recommended when an adequate retention is suitable, and the preparations have big parallel walls.

CLINICAL BENEFITS:

- Avoid allergic reactions or sensitivity.
- Easy removal.
- Compatibility with modern restorations.
- Good temporary sealing.
- Rapid initial hardening (60-80 sec.)
- It does not stain temporary acrylic teeth.

PERFORMANCE CHARACTERISTICS:

- The permanent cementation doesn't have to exceed the estimated time to use non-eugenol cement.
- They have lower retention than other cements or resin cements, especially when Retrievability is desired.
- Mechanical and chemo-mechanical cleaning protocols of non-eugenol cements can influence restoration bond strength. Clinically relevant sealing for provisional restorations in short-term simulations.
- Microleakage in Zinc Oxide Non-eugenol is more in comparison with others cements. Adequate mechanical resistance. More homogeneous texture. Reduced solubility. Moderate short-term compressive strength.

COLOR: Dentin color (light beige/cream).

INTENDED USE: Provisional cementation of temporary dental restorations or prostheses, such as crowns, bridges, or inlays, during the period prior to the definitive cementation.

INDICATIONS:

- When planning to permanently cement with adhesive resins.
- In patients with high aesthetic demands or moderate masticatory pressure.
- As a temporary cement in cases of adhesive indirect restorations such as inlays, veneers and crowns.
- When it is necessary to fully preserve the adhesive potential of dentin and the internal restoration substrate.
- In orthodontic cases where temporary crowns are required for several months.
- When occlusion needs adjusting.
- In patients with sensitivity or allergic reactions to eugenol.

INTENDED USER: Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must be make restorations with Medental Temporary Cement (Eugenol Free).

PRECAUTIONS: Some patients may present allergy to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

Does not require refrigeration.

Only for professional dental use.

Room temperature before use. It does not require cold chain to transport.

Do not use after date of expiry.

CONTRAINDICATIONS:

- Do not use if the tooth is very sensitive post-operatively.
- If a sedative effect is sought in a preparation that causes discomfort.
- If it is required that the cement be removed more easily.
- If long-term permanent cementation is required.
- If pulp irritation exists in deep cavities.
- Postoperative discomfort in deep areas.
- Avoid in crowns, bridges, posts and definitive or permanent restorations.
- Avoid use in highly sensitive pediatric patients.
- If there is an allergy to any of the components of the formula.

COMPOSITION: *Base:* Zinc oxide, Peg and Stabilizers. *Catalyst:* Resin, highly dispersed Silica Dioxide, Isobutyl Salicylate, Plastifying, Stabilizers and Catalyst.

INSTRUCTIONS FOR USE:

A. GENERAL INFORMATION: Temporary Cement (Eugenol Free) has a base made out of Zinc Oxide and the Catalyst has a resin which reacts with the Zinc Oxide to create cement with a strong bond. The cementation layer is thin enough to avoid interfering with the suitable setting of the prosthesis.

B. OUTSTANDING FEATURES: - Non-eugenol zinc oxide-based cement. - No irritation of the pulp and has good compatibility. - Initial hardness fast (60-80 sec) - Does not inhibit the hardening reaction of permanent cements or restorative materials made out of resins. - Bonds to major temporary materials containing methyl-methacrylate. -Does not stain temporary acrylics. - Cement without odor, dentin color.

C. PHYSICAL PROPERTIES: Compressive Strength: 35MPa (after 24 hour) Film thickness: 15 microns Working time: 60-80 sec. max. Setting time: 5 min.

D. INSTRUCTIONS FOR CEMENTATION:

1. Clean the tooth preparation with water; be sure there is no residual dentine.
2. In areas close to the pulp, must apply a thin layer of pulp protector.
3. Damaged areas of the preparation of the tooth can be rebuilt with glass ionomer, compomer or composite.

4. Put mixing tip on the double syringe, lock it and dispose the cement directly onto the prosthesis, after set firmly the prosthesis upon the tooth or teeth.

5. Be sure the cement has a "beige" color and not white. This is indicative that base and catalyst are well mixed.

6. Be sure there is no residual cement caught under the gingival margins or under the prosthesis.

7. Working time is less than 80 sec. depending on the heat and humidity relatives to the practical location, total hardening will occur in 5 min. after the product has been exposed to the body temperature.

8. Temporary cement (Eugenol Free) can be removed from the gingival margins after 80 sec. corresponding to the working time.

CONSERVATION:

- Close the double syringe immediately after dispensing the product.
- Always store the product in a fresh area.

WARNINGS and PRECAUTIONS: In case of any change on the original function of the product, stop using it and contact the manufacturer.

- Read carefully and follow all instructions.

- It may cause mild pulp irritation or inflammation immediately after its application.

- It can cause local allergies and contact dermatitis. Its clinical use is relatively safe if handled properly.

- It is recommended to thoroughly clean the remnants before using the final cement or resin.

- Should not be used as a permanent cement.

- Should not be used if a strong bond is required.

- In deep cavities, place a protective intermediate base to avoid irritation.

- Before the use of these products always use rubber dam gloves and protective glasses.

- In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor.

- In case of contact with skin wash the area with soap and water.

- The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately

and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance.

- Ask for patient to rinse his mouth immediately after treatment.

- The product is designed to be used at room temperature.

WARRANTY: If the use of this product on a patient or anyone else causes some anomaly, please send detailed report to: Conamco, SA de CV email:

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Warranty: MEDENTAL INTERNATIONAL and/or CONAMCO SA de CV warranty the product against any fault of manufacturing of the material. MEDENTAL INTERNATIONAL AND/OR CONAMCO SA de CV does not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product. The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at: <https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS: New, clean nitrile gloves should be worn.

1) *Removal of excess.*

- Fresh (before complete setting): Remove with a probe or spatula. recommends cleaning immediately.

- Setting: Using cement chisels, curettes or metal excavators with controlled movements, pass waxed dental floss in interproximal areas, pulling downwards to avoid fracturing restorations.

2) *Removal under restorations or crowns.*

- It is not possible to remove it without removing the restoration; in these cases, the crown must be sectioned with a carbide/diamond bur and the tooth cleaned manually or with ultrasound.

*In both cases, it is recommended to protect soft tissues, work with good lighting, and use a magnifying glass to identify cement residues.

3) Dispose of with solid waste in accordance with the regulations of each country.

FEATURES CLAIMS: Zinc Oxide Free-Eugenol based cements are used as temporary cements short-term. Eugenol-free temporary cements are effective in the clinical setting in the short term. Since they are not used as permanent solutions, they are easy to remove, but their effectiveness is short-term, so it is recommended not to leave too much time between applications. Temporary and permanent restoration. It is necessary to ensure that no remnants remain. to ensure the adhesion of the final restorative material. To take good care of the restoration, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and use dental protector if the patient has Bruxism.

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