MULTIPURPOSE ZOE CEMENT

INSTRUCTIONS FOR USE

English
Español
Français

 CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd.-1st.
P.C. 08025, Barcelona, Spain.

REV: "0" / 21 JUL 2025

MULTIPURPOSE ZOE CEMENT

ENGLISH

PRODUCT DESCRIPTION:

ZOE Cements are a family of products made with Zinc oxide powder (some with strengtheners and accelerators) and eugenol liquid. Eugenol is the essence of cloves. ZOE cements are generally used as temporary cements with expected duration. Multipurpose ZOE cement is a dressing/surgical cement to serve as a protective covering over operated or treated areas to prevent contact with irritating or infecting substances.

CLINICAL BENEFITS:

- Eugenol has a sedative effect that relieves pain and inflammation.
- Eugenol inhibits the formation of bacteria.
- It has good thermal insulation.
- It is economical and easy to use.

PERFORMANCE CHARACTERISTICS:

- Good sealing in temporary restorations, helping to protect the dentin/pulp.
- Its mechanical properties (low resistance) facilitate removal.
- It has no negative effect on bond strength (only if used with light-curing resin).

COLOR: Transparent

INTENDED USE: For temporary and intermediate cementation and filling.

INDICATIONS:

- Dressing/Surgical to serve as a protective covering over operated or treated areas to prevent contact with irritating or infecting substances.
- Temporary filling material for cavities.
- Temporary luting Cement for the cementation of temporary crowns and bridges.

INTENDED USER: Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must make restorations with Medental Multipurpose ZOE Cement.

PRECAUTIONS: Some patients may present allergy to the components of the formula, in this case, suspend its use.

Keep away from Sunlight.

Keep away from rain.

Keep away the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

Does not require refrigeration.

Only for professional dental use.

Room temperature before use. It does not require cold chain to transport.

Do not use after date of expiry.

CONTRAINDICATIONS:

- The material described must not be used when the application techniques are not possible.
- It is recommended not to apply or remove it completely if light-curing resins are used, since eugenol inhibits polymerization.
- It is not recommended if a strong and long-lasting bond is required.
- Do not use in people allergic to eugenol.
- Do not use in deep cavities near the nerve.

COMPOSITION: Zinc Oxide, Talc, Resinall, Eugenol, Propanol resinous oil and camphor synthetic.

INSTRUCTIONS FOR USE:

- A. SURGICAL DRESSING:** The powder and liquid should be mixed until a non-sticky homogeneous paste is obtained. The paste is then applied and should be left in place in sockets, on dentures, over treated or operated areas, a required number

of days as suggested in the techniques. This will allow sufficient time for healing. Premature removal will disturb new granulations. However, there is a difference of opinion among dentists as to the length of time temporary cement should be left in place. However, in general it should be replaced with a new one after 7-10 days. When the use of Temporary cement is recommended in sockets or as a covering over wounds it is understood that the material must be removed by the dentist. If a socket, contrary to the directions, were packed tightly with gauze or cotton impregnated with Multi-purpose ZOE Cement and the packing left in place too long, healing would materially be delayed. Consequently, the packing should not be left in place longer than specified time, maximum two weeks, and the patients should always be instructed to return to have the pack removed. Multi-purpose ZOE Cement is manufactured to set slowly out of the mouth, so as to give sufficient time for manipulation as is required in periodontal surgery for instance, but sets quickly in the mouth on contact with blood or saliva.

- B. MIXING:** Whether mixed thin or applied on gauze into sockets or mixed thick as required following periodontal surgery. Multi-Purpose Cement should be thoroughly mixed. The powder must be properly incorporated into the liquid and thoroughly spatulated. Take time and mix properly. The pack should also be thoroughly spatulated to ensure uniform setting. The proportions indicated in the different techniques shown below are approximations. The accurate quantity of powder and liquid depends on the consistency desirable by the dentist.
- C. EXTRACTIONS:** Mix thin (0.8g of powder by 0.5ml of liquid) Multi-purpose ZOE Cement and apply on a strip of gauze or cotton string and place lightly into the socket so as to form a support for

the blood clot. Do not pack tightly. Leave in place 3 to 7 days.

- D. DRY SOCKETS:** Carefully remove debris and fill socket with Multi-Purpose ZOE Cement Paste. Mix cement thin (0.8g of powder with 0.5ml of liquid) and apply on a strip of gauze or cotton string and place lightly into the socket. Leave in place 3 to 7 days if patient is comfortable. It may be necessary to replace the packing. In which case repeat the procedure described, but it is not advisable to renew the pack for at least 24 hours. Do not pack lightly.
- E. PERIODONTAL SURGERY:** Following surgery mix Multi-purpose Cement very thick (2.2g of the powder with 5ml of liquid) and apply interproximally, packing toward the contact points and force the pack through onto the lingual surfaces. Cover all operated areas and lock in place by pressure of the thumb and forefinger. Leave in place 10 to 14 days.
- F. VICENT'S INFECTION:** Treat and medicate as is your custom thin pack over treated areas with thick mix (2.2g of powder with 0.5ml of liquid) of Multi-purpose ZOE Cement into which has been incorporated 1/6 part of silver iodide powder. Pack as, directed for periodontal surgery interproximally. If possible, leave in place 7 to 10 days.
- G. TEMPORARY FILLINGS:** Mix (2g of powder with 0.5ml of liquid) of Multi-purpose ZOE Cement until a thick firm paste is obtained. Incorporate a little cotton if desirable and pack tightly into the cavity.
- H. IMMEDIATE DENTURES:** Mix (2g of powder with 0.5ml of liquid) Multi-Purpose ZOE Cement until a heavy, but soft enough that it will flow under pressure. Instruct patient in usual cheek, lip and tongue manipulations and maintain bite firmly until pack sets. Add a few drops of water to accelerate setting.

I. TEMPORARY REBASES: Mix and apply as for immediate dentures.

J. TEMPORARY LUTING: Mix (0.8g of powder by 0.5ml of liquid) of Multi-purpose ZOE Cement until a thin creamy paste is obtained and apply freely into the cavity and on the inside surfaces of the temporary inlay, bridge or crown. Before applying Multi-purpose into temporary bridges or inlays, place a thin film of wax on the cavity surface of the bridge, crown or inlay. It can then be easily removed by first heating with a hot instrument to soften the wax lining.

CONSERVATION:

- Close the package immediately after dispensing the product.
- To prevent deterioration of material, keep bottle and jar tightly closed when not in use.

WARNINGS and PRECAUTIONS: In case of any change on the original function of the product, stop using it and contact the manufacturer.

- Read carefully and follow all instructions.
- It may cause mild pulp irritation or inflammation immediately after its application.
- It can cause local allergies and contact dermatitis. Its clinical use is relatively safe if handled properly.
- It is recommended to thoroughly clean the remnants before using the final cement or resin.
- Should not be used as a permanent cement.
- Should not be used if a strong bond is required.
- In deep cavities, place a protective intermediate base to avoid irritation.
- Before the use of these products always use rubber dam gloves and protective glasses.
- In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor.
- In case of contact with skin wash the area with soap and water.
- The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately

and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance.

- Ask for patient to rinse his mouth immediately after treatment.

- The product is designed to be used at room temperature.

WARRANTY: If the use of this product on a patient or anyone else causes some anomaly, please send detailed report to: Conamco, SA de CV email:

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Warranty: MEDENTAL INTERNATIONAL and/or CONAMCO SA de CV warranty the product against any fault of manufacturing of the material. MEDENTAL INTERNATIONAL AND/OR CONAMCO SA de CV does not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product. The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at:

<https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS:

New, clean nitrile gloves should be worn.

1) *Removal of excess.*

- Fresh (before complete setting): Remove with a probe or spatula. recommends cleaning immediately.
- Setting: Using cement chisels, cures, or metal excavators with controlled movements, pass waxed dental floss in interproximal areas, pulling downwards to avoid fracturing restorations.

2) *Removal under restorations or crowns.*

- It is not possible to remove it without removing the restoration; in these cases, the crown must be sectioned with a carbide/diamond bur and the tooth cleaned manually or with ultrasound.

*In both cases, it is recommended to protect soft tissues, work with good lighting, and use a magnifying glass to identify cement residues.

3) Dispose of with solid waste in accordance with the regulations of each country.

FEATURES CLAIMS: Zinc oxide and eugenol-based cements are used as short- and medium-term temporary cements. ZOE cements are effective in the short- and medium-term clinical setting (reinforced with hardener). Their main characteristic is their sedative and antibacterial effect. Since they are not used as permanent cements, they are easy to remove, but their performance is short-term, so it is recommended not to leave too much time between the temporary and permanent restorations. It is necessary to ensure that no excess fluid remains to ensure adhesion of the final restorative material, especially if they are light-curing resins. To take good care of the restoration, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and dental use protector if the patient has Bruxism.

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