## Z.O.E. ROOT CANAL SEALER

### INSTRUCTIONS FOR USE

English  
Español  
Français

 CONAMCO, S.A. de C.V.  
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.  
Col. Algarín, Alcaldía Cuauhtémoc.  
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.  
Travessera de Gracia Street No. 315, 2nd.-1st.  
P.C. 08025, Barcelona, Spain.

Rev:"0" / 21 JUL 2025

## Z.O.E. ROOT CANAL SEALER

### ENGLISH

**PRODUCT DESCRIPTION:** ZOE Root Canal Sealer cement is a zinc oxide eugenol cement with a powder reinforcement and barium sulfate to make it radiopaque. It is a cement to fill the root canal in periapical surgery.

**CLINICAL BENEFITS:**

- Eugenol has a sedative effect that relieves pain and inflammation.
- Eugenol inhibits the formation of bacteria.
- It has good thermal insulation.
- It is economical and easy to use.

**PERFORMANCE CHARACTERISTICS:**

- Good sealing in temporary restorations, helping to protect the dentin/pulp.
- Its mechanical properties (low resistance) facilitate removal.
- It has no negative effect on bond strength (only if used with light-curing resin).

**COLOR:** Transparent

**INTENDED USE:** For temporary and intermediate cementation and filling.

**INDICATIONS:**

For the cementation of the gutapercha during endo treatment.

**INTENDED USER:** Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

**PATIENT POPULATION:** All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must make restorations with Medental ZOE Root Canal Cement.

**PRECAUTIONS:** Some patients may present allergy to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

Does not require refrigeration.

Only for professional dental use.

Room temperature before use. It does not require cold chain to transport.

Do not use after date of expiry.

**CONTRAINDICATIONS:**

- The material described must not be used when the application techniques are not possible.
- It is recommended not to apply or remove it completely if light-curing resins are used, since eugenol inhibits polymerization.
- It is not recommended if a strong and long-lasting bond is required.
- Do not use in people allergic to eugenol.
- Do not use in deep cavities near the nerve.

**COMPOSITION:** Zinc Oxide, Eugenol, Bismuth Subnitrate, Rosin lump.

**INSTRUCTIONS FOR USE:**

**Note:** The ZOE Root Canal Sealer Cement is radiopaque for the purpose of sealing of root canal.

A. Place two (2) drops of the liquid on a sterile slat for two (2) spoons of powder.

B. For precise dosage and lab test, weight 0.5g of powder for 0.1ml. of liquid.

C. Spatula one portion into the liquid until all particles of powder are dissolved. Add more portions of powder to the mix, spatulating until a creamy paste, very smooth are obtained. Spatulating time ranges from 3 to 5 minutes. The longer spatulating time, the faster the setting time. When the mixed cement is of such consistency that a thin unbroken thread of at least 1 inch lengths spatulation is completed, may be obtained when the flat surface of the spatula is slowly raised from the cement and slab.

D. Of the root canal cement, the final setting time on the glass slab a relative humidity of approximately 40% at 72°F room temperature should be between 8 ½ and 10 hours. In the root canal such a cement mixture will have an initial setting time approximately 20 minutes, thus enabling the Dentist to make a correction should it be necessary.

**CONSERVATION:**

- Close the package immediately after dispensing the product.
- To prevent deterioration of material, keep bottle and jar tightly closed when not in use.

**WARNINGS and PRECAUTIONS:** In case of any change on the original function of the product, stop using it and contact the manufacturer.

- Read carefully and follow all instructions.
- It may cause mild pulp irritation or inflammation immediately after its application.
- It can cause local allergies and contact dermatitis. Its clinical use is relatively safe if handled properly.
- It is recommended to thoroughly clean the remnants before using the final cement or resin.
- Should not be used as a permanent cement.
- Should not be used if a strong bond is required.
- In deep cavities, place a protective intermediate base to avoid irritation.
- Before the use of these products always use rubber dam gloves and protective glasses.
- In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor.
- In case of contact with skin wash the area with soap and water.
- The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance.
- Ask for patient to rinse his mouth immediately after treatment.
- The product is designed to be used at room temperature.

**WARRANTY:** If the use of this product on a patient or anyone else causes some anomaly, please send detailed report to: Conamco, SA de CV email:

[atencionaclientes@conamcosacv.com.mx](mailto:atencionaclientes@conamcosacv.com.mx)

Warranty: MEDENTAL INTERNATIONAL and/or CONAMCO SA de CV warranty the product against any fault of manufacturing of the material. MEDENTAL INTERNATIONAL AND/OR CONAMCO SA de CV does not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product. The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at:

<https://ec.europa.eu/tools/eudamed>

**DISPOSAL CONSIDERATIONS:**

New, clean nitrile gloves should be worn.

1) *Removal of excess.*

- Fresh (before complete setting): Remove with a probe or spatula. recommends cleaning immediately.
- Setting: Using cement chisels, curettes or metal excavators with controlled movements, pass waxed dental floss in interproximal areas, pulling downwards to avoid fracturing restorations.

2) *Removal under restorations or crowns.*

- It is not possible to remove it without removing the restoration; in these cases, the crown must be sectioned with a carbide/diamond bur and the tooth cleaned manually or with ultrasound.

\*In both cases, it is recommended to protect soft tissues, work with good lighting, and use a magnifying glass to identify cement residues.

3) Dispose of with solid waste in accordance with the regulations of each country.

**FEATURES CLAIMS:** Zinc oxide and eugenol-based cements are used as short- and medium-term temporary cements. ZOE cements are effective in the short- and medium-term clinical setting (reinforced with hardener). Their main characteristic is their

sedative and antibacterial effect. Since they are not used as permanent cements, they are easy to remove, but their performance is short-term, so it is recommended not to leave too much time between the temporary and permanent restorations. It is necessary to ensure that no excess fluid remanent to ensure adhesion of the final restorative material, especially if they are light-curing resins. To take good care of the restoration, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and dental use protector if the patient has Bruxism.

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