## Z.O.E. CEMENT TYPE III

### INSTRUCTIONS FOR USE

English  
Español  
Français

 CONAMCO, S.A. de C.V.  
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.  
Col. Algarín, Alcaldía Cuauhtémoc.  
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.  
Travessera de Gracia Street No. 315, 2nd.-1st.  
P.C. 08025, Barcelona, Spain.

Rev: "0" / 21 JUL 2025

## Z.O.E. CEMENT TYPE III

### ENGLISH

**PRODUCT DESCRIPTION:** ZOE Cements are a family of products made with Zinc oxide powder (some with strengtheners and accelerators) and eugenol liquid. Eugenol is the essence of cloves. ZOE cements are generally used as temporary cements with expected duration. The ZOE Type III cement is the longest lasting ZOE temporary cement. Es reinforced with acrylic powder.

#### CLINICAL BENEFITS:

- Eugenol has a sedative effect that relieves pain and inflammation.
- Eugenol inhibits the formation of bacteria.
- It has good thermal insulation.
- It is economical and easy to use.

#### PERFORMANCE CHARACTERISTICS:

- Good sealing in temporary restorations, helping to protect the dentin/pulp.
- Its mechanical properties (low resistance) facilitate removal.
- It has no negative effect on bond strength (only if used with light-curing resin).

**COLOR:** Transparent

**INTENDED USE:** For temporary and intermediate cementation and filling.

**INDICATIONS:** For odontological intermediate restorations particularly odontopediatrics. For use with rampant caries where the demineralized tissue should all be removed from the cavity as fast as possible to reduce the concentration of cariogenic bacteria and prevent the process that produces caries, during this time the teeth ought to be protected by ZOE cement type III which is a restoration of intermediate duration. It can also be used for bases or as a temporary cement,

it also has the advantage of protection and healing the pulp.

**INTENDED USER:** Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

**PATIENT POPULATION:** All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must make restorations with Medental ZOE Type III Cement.

**PRECAUTIONS:** Some patients may present allergy to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

Does not require refrigeration.

Only for professional dental use.

Room temperature before use. It does not require cold chain to transport.

Do not use after date of expiry.

**CONTRAINDICATIONS:** • The material described must not be used when the application techniques are not possible. • It is recommended not to apply or remove it completely if light-curing resins are used, since eugenol inhibits polymerization. • It is not recommended if a strong and long-lasting bond is required. • Do not use in people allergic to eugenol. • Do not use in deep cavities near the nerve.

**COMPOSITION:** Zinc Oxide, Eugenol, Hardener, Acrylic Polymer.

#### INSTRUCTIONS FOR USE:

For odontological intermediate restorations particularly odontopediatrics. For use with rampant caries where the demineralized tissue should all be removed from the cavity as fast as possible to reduce the concentration of cariogenic bacteria and prevent the

process that produces caries, during this time the teeth ought to be protected by ZOE Cement Type III which is a restoration of intermediate duration. It can also be used for bases or as temporary cement; it also has the advantage of protecting and healing the pulp. Note: ZOE Cement Type III can be used for cementation of intermediate restorations (two to three months) of crowns and bridges and temporary not definitely inlays. Before The use of this product always use rubber dam to avoid any spill or splash over the teeth.

A. Take a portion of powder in leveled scoops sufficient to produce, with the liquid, a mix of smooth and flexible consistency to fill the cavity.

B. Place the powder on a slab and divide into two parts. Next divide these two parts into a further two parts each.

C. On another part of the slab place 1 drop of liquid each leveled scoop.

D. With a stainless-steel spatula, mix each portion of powder with the liquid thoroughly.

E. Mixing each part for 20-30 seconds. Total time of mixing about 80 seconds.

F. Place the mixture in the previously prepared cavity, dried and with mechanical retention. Apply pressure to permit internal adjustment by asking the patient to bite firmly for mix is smooth before adding more powder.

G. The mix will be complete and rather stiff in 1 ½ minutes.

#### **CONSERVATION:**

- Close the package immediately after dispensing the product.

- To prevent deterioration of material, keep bottle and jar tightly closed when not in use.

**WARNINGS and PRECAUTIONS:** In case of any change on the original function of the product, stop using it and contact the manufacturer.

- Read carefully and follow all instructions.

- It may cause mild pulp irritation or inflammation immediately after its application.

- It can cause local allergies and contact dermatitis. Its clinical use is relatively safe if handled properly.

- It is recommended to thoroughly clean the remnants before using the final cement or resin.

- Should not be used as a permanent cement.

- Should not be used if a strong bond is required.

- In deep cavities, place a protective intermediate base to avoid irritation.

- Before the use of these products always use rubber dam gloves and protective glasses.

- In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor.

- In case of contact with skin wash the area with soap and water.

- The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance.

- Ask for patient to rinse his mouth immediately after treatment.

- The product is designed to be used at room temperature.

**WARRANTY:** If the use of this product on a patient or anyone else causes some anomaly, please send detailed report to: Conamco, SA de CV email:

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Warranty: MEDENTAL INTERNATIONAL and/or CONAMCO SA de CV warranty the product against any fault of manufacturing of the material. MEDENTAL INTERNATIONAL AND/OR CONAMCO SA de CV does not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product. The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at:

<https://ec.europa.eu/tools/eudamed>

#### **DISPOSAL CONSIDERATIONS:**

New, clean nitrile gloves should be worn.

1) *Removal of excess.*

- Fresh (before complete setting): Remove with a probe or spatula. recommends cleaning immediately.

- Setting: Using cement chisels, cures or metal excavators with controlled movements, pass waxed dental floss in interproximal areas, pulling downwards to avoid fracturing restorations.

2) *Removal under restorations or crowns.*

- It is not possible to remove it without removing the restoration; in these cases, the crown must be sectioned with a carbide/diamond bur and the tooth cleaned manually or with ultrasound.

\*In both cases, it is recommended to protect soft tissues, work with good lighting, and use a magnifying glass to identify cement residues.

3) Dispose of with solid waste in accordance with the regulations of each country.

**FEATURES CLAIMS:** Zinc oxide and eugenol-based cements are used as short- and medium-term temporary cements. ZOE cements are effective in the short- and medium-term clinical setting (reinforced with hardener). Their main characteristic is their sedative and antibacterial effect. Since they are not used as permanent cements, they are easy to remove, but their performance is short-term, so it is recommended not to leave too much time between the temporary and permanent restorations. It is necessary to ensure that no excess fluid remanent to ensure adhesion of the final restorative material, especially if they are light-curing resins. To take good care of the restoration, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and dental use protector if the patient has Bruxism.

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