POLYCARBOXYLATE CEMENT

P.C.A.

INSTRUCTIONS FOR USE

English
Español
Français

CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd.-1st.
P.C. 08025, Barcelona, Spain.

Rev: "0" / 21 Jul 2025

POLYCARBOXYLATE CEMENT (PCA)

ENGLISH

PRODUCT DESCRIPTION:

Polycarboxylate Cement is a material used for cementations of inlays, crowns, bridges and orthodontic appliances.

CLINICAL BENEFITS:

- Long history of use
- Excellent compressive strength making the restoration more durable
- Allows for thin film thickness making it more adaptable and precise.
- Once set, it remains stable.
- Chemical adhesion to enamel and dentin.
- Less irritant than Zinc Phosphate Cement.
- Working time.
- Compatible with metal and ceramic.
- Antibacterial effect.
- It returns the masticatory forces of the treated piece and in cases or orthodontic keeps the devices in position during the treatment.
- Rebuild and return chewing functionality.

PERFORMANCE CHARACTERISTICS:

- High tensile strength.
- Chemical Adhesion.
- Moderate solubility.
- Moderate working time.
- More Biocompatible.
- Less mechanical resistance.
- More solubility than other cements.
- Not indicated as temporary cement.

COLOR: Opaque.

INTENDED USE:

For permanent luting cementation of crowns and bridges restorations. (To attach or cement indirect restorations to prepared teeth).

INDICATIONS:

- Permanent cementation of crowns, inlays, bridges and orthodontic brackets.
- Cavitary base beneath amalgam and metal restorations, especially in teeth with a large pulp chamber (less irritating than Zinc Phosphate)
- Cavitary lining in deep cavities.
- As an alternative to Zinc Phosphate when a chemical bond to the enamel and dentin.

INTENDED USER: Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must make restorations with Medental Polycarboxylate Cement.

PRECAUTIONS:

Some patients can present allergies to the components of the formula, in this case, suspend its use.

Keep away from Sunlight.

Keep away from rain.

Keep away from the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

No requires refrigeration.

Only for professional dental use.

Room temperature before use.

Not require cold chains to transport.

Do not use it after the date of expiry.

CONTRAINDICATIONS: • The material described must not be used when the application techniques are not possible. • It is not indicated for use as a temporary cement because it can be difficult to removal. • Very large cavities or those with thin walls due to their low resistance. • It is opaque, so it is not recommended when there is a need for great aesthetics. • Cases where greater fluidity or manipulation is required, since their working time it's short. • Patients with a history of

allergic reactions to polyacrylic acid • Areas that require high mechanical resistance.

COMPOSITION: Zinc Oxide and Polyacrylic acid.

INSTRUCTIONS FOR USE:

NOTE: If the cavity is deep enough to compromise the pulp is recommended to use a layer of calcium hydroxide before the use of this cement.

A. Dispensing Material:

1. Take powder from the jar with the dispensing spoon, is the quantity required and put it on a mixing pad or glass slab.

2. Take the liquid bottle and dispense the number of drops required.

B. Mix Ratio: For luting depending on the consistency required, the quantity for cementation recommended is: One spoon and one drop, but the dentist can change to others proportions. Or lining and filling, the quantity recommended are 3 spoons of powder and 2 drops of liquid, but the dentist can change to other proportions. To test the product in a lab, the quantity needed is: 1.0 g of powder and 0.5 g of liquid.

C. Mixing: Using a suitable incorporate all the powder rapidly into the liquid. Mix until a creamy consistency is reached for a cementation material, or a putty-like consistency in case of a lining or filling needs.

D. Mixing Time: 30 seconds. **Working Time:** 2 minutes (from end of mix).

F. Application: Apply cement immediately on the completion of the mix until the appearance is glossy. Do not use when shiny aspect is lost.

G. Finishing: Excess material may be removed at once or when the cement has completely hardened. Do not touch the material when it is in the intermediate rubbery condition.

H. Setting Time: 6-7 minutes.

CAUTION: Due to the adhesive nature of Polycarboxylate Cement, spatulas, glass slabs and placement instruments should be wiped clean under running water immediately after use.

CONSERVATION:

- Close the package immediately after dispensing the product. - To prevent deterioration of material, keep bottle and jar tightly closed when not in use. - To eliminate condensation of water vapor on the slab, raise the temperature of the slab slightly if 72°F (22°C) temperature of slab is below the dew point. - Be sure all the powder is thoroughly incorporated into the liquid in each successive step of the mixing.

WARNINGS and PRECAUTIONS: In case of any change to the original function of the product, stop using it and contact the manufacturer.

- Read carefully and follow all instructions - May cause mild pulp irritation or inflammation during application. - It can be difficult to remove if used as temporary cement. - Before the use of these products always use rubber dam gloves and protective glasses. - In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor. - In case of contact with skin wash the area with soap and water. - The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance. - Ask for patient to rinse his mouth immediately after treatment. - The product is designed to be used at room temperature.

WARRANTY:

If the use of this product on a patient or anyone else causes some anomaly, please send detailed report to: Conamco, SA de CV email:

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Warranty: MEDENTAL INTERNATIONAL and/or CONAMCO SA de CV warranty the product against any fault of manufacturing of the material. MEDENTAL INTERNATIONAL AND/OR CONAMCO SA de CV does not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product. The summary of Safety and clinical

Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at: <https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS:

New, clean nitrile gloves should be worn.

1) Removal of excess.

- Fresh (before complete setting): Remove with a probe or spatula. recommends cleaning immediately.

- Setting: Using cement chisels, curettes or metal excavators with controlled movements, pass waxed dental floss in interproximal areas, pulling downwards to avoid fracturing restorations.

2) Removal under restorations or crowns.

- It is not possible to remove it without removing the restoration; in these cases, the crown must be sectioned with a carbide/diamond bur and the tooth cleaned manually or with ultrasound.

* In both cases, it is recommended to protect soft tissues, work with good lighting, and use a magnifying glass to identify cement residues.

3) Dispose of with solid waste in accordance with the regulations of each country.

FEATURES/CLAIMS: Water Base Cements are a material used for cementations of inlays, crowns, bridges and orthodontic appliances and occasionally as a temporary restoration. It is prepared by mixing Zinc Oxide and magnesium oxide powders with a liquid consisting mainly of polyacrylic acid, water and buffers. This Cements can be used for fixing crowns, bridges, posts, cores, in various materials with clinical properties like modulus resilience, mechanical adhesive, solubility, adhesion and retention, long term survival event rates. Water Base cements are effective in clinical context and in short-term clinical use, and they can increase strength with storage. Polycarboxylate Cement in applications emphasizing bonding and seal integrity, it is more flexible, resilient, adhesive and has a good seal for microleakage. To take good care of the restoration, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as

coffee and wine, as well as maintain good dental hygiene and use dental protector if the patient has Bruxism.

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