ETCHANT

(37% ORTHOPHOSPHORIC ACID)

INSTRUCTIONS FOR USE

English
Español
Français

 CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd.-1st.
P.C. 08025, Barcelona, Spain.

Rev: "0" / 21 JUL 2025

ETCHANT (37% ORTHOPHOSPHORIC ACID)

ENGLISH

DESCRIPTION: The MEDENTAL Etchants are designed as part of a system to enhance the bonding of composite resin for tooth restoration. They are formulated with orthophosphoric acid at a concentration of 37% and are available in gel or liquid form for easy application. Orthophosphoric acid produces a controlled surface demineralization in the tooth enamel, making it rough to facilitate the penetration of adhesive systems.

CLINICAL BENEFITS:

- Strong chemical adhesion
- Long-term duration.
- Easy to use
- Stays in place and will not run onto exposed tissues
- Washes off easily with water

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SHADES: Blue or green gel

INTENDED USE: For conditioning enamel and dentin before the application of adhesives and restorative materials.

INDICATIONS: Etchants can be used for a wide range of clinical applications today, because they are part of many dental treatments with composites, resin cements, sealants, etc. Etchants are indicated:

- Repair of resin restorations.
- Placement of orthodontic brackets.
- Cementation of indirect restorations.
- Application of the pit and fissures sealant

INTENDED USER: Dental professionals only with an academic degree, experience, and training for use

dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must make restorations with Medental Etchants.

PRECAUTIONS:

Some patients can present allergies to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away from the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

No requires refrigeration.

Only for professional dental use.

Room temperature before use.

Not require cold chains to transport.

Do not use it after the date of expiry.

CONTRAINDICATIONS:

- The material must not be used when the application techniques described are not possible.
- Hypersensitivity to the components of the formula.
- Do not use if the dentin is close to the pulp or if there is pulp exposure, as it may cause irritation.
- Do not use on hypersensitive lesions
- Patients with multiple allergies
- Patients with ulcerative lesions in the mouth
- Patients with respiratory diseases

COMPOSITION: Phosphoric acid, bidistilled water, pigment and gelling agent.

INSTRUCTIONS FOR USE:

1. Preparing the teeth

Once identified, the tooth is properly isolated, using staples and a medium-thickness rubber dam.

2. Enamel engraving

Enamel and dentin are etched with 35% orthophosphoric acid for 20 to 30 seconds on enamel

and 10 to 15 seconds on dentin. It is cleaned with water and it is confirmed that the enamel already has a greyish appearance, not drying out the tissue and keeping it moist. Otherwise, the operation is repeated.

Clinical Procedure:

- a) Clean and dry the cavity and tooth.
- b) Etch the enamel of the tooth with the etchant acid for 20 or 30 sec.
- c) Rinse and dry with clean air.
- d) If the enamel does not have a rough chalk aspect repeat the process.

3. Continue with the desired procedure:

- Adhesive, Resin, Sealant or other.

- Once the etching time is over, rinse with plenty of water for at least 10-15 seconds to completely remove the acid. - Aspirate excess water with the diluted acid and avoid contamination of soft tissues. - Dry if it is enamel (if it is dentin it should remain slightly damp) - Assess whether irritation persists or if sensitivity is present.

CONSERVATION:

- The system is designed to be used at room temperature above 24°C or 75°F.
- If stored in refrigerator, bring to room temperature prior to use.
- Do not expose materials to elevated temperature.

WARNINGS and PRECAUTIONS:

-Read carefully and follow all instructions
-Before the use of these products always use rubber dam gloves and protective glasses.
-The etching solution contains a 37% phosphoric acid solution, which can be harmful to skin and eyes.
-In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor.
-In case of contact with skin wash the area with soap and water.
-The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately

and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance.

-Ask for patient to rinse his mouth immediately after treatment.

-The product is designed to be used at room temperature. If refrigerated, the material may be difficult to extrude.

- Do not exceed the engraving time.

WARRANTY:

If the use of this product presents some anomaly such as:

- It does not fulfill the intended purpose,
- It doesn't return the functionality of the tooth,
- It doesn't have the properties it refers to,
- It has damage to packaging or manufacturing,
- The instructions for use are misunderstood,
- Any secondary reaction or anomaly:

Please send detailed report to: Conamco, S.A. de C.V. email: atencionclientes@conamcosacv.com.mx

Warranty: MEDENTAL and/or CONAMCO S.A. de C.V. warranty the product against any fault of manufacture of the material. MEDENTAL and/or CONAMCO S.A. de C.V. do not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective products.

The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at:

<https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS:

-New, clean nitrile gloves should be worn.
- To discard it: It is necessary to dilute it with water and not throw it down the drain.

FEATURES/CLAIMS:

- Improves resin infiltration.
- It is a biocompatible material if used correctly
- High mechanical resistance.
- Excellent adhesion.
- Engraving depth 5-50 µm.

- Clinically versatile.
- It is the reference standard for enamel adhesion.
- The best agent for a most effective enamel etching.

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-Veuillez envoyer un rapport détaillé à : Conamco, S.A. de C.V., email :