DUAL CURE RESIN CEMENT

INSTRUCTIONS FOR USE

English
Español
Français

 CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd. -1st.
P.C. 08025, Barcelona, Spain.

DUAL CURE RESIN CEMENT

ENGLISH

PRODUCT DESCRIPTION

DESCRIPTION: Dual-cure resin is a dental material designed to polymerize by both light and chemical reaction, ensuring polymerization where light cannot fully penetrate.

CLINICAL BENEFITS:

- Working time
- Strong chemical adhesion
- Long-term duration.
- Easy to use (fewer clinical steps)

* Assured Polymerization

PERFORMANCE CHARACTERISTICS:

- High mechanical resistance
- Excellent adhesion.
- Long-term stability.

SHADES: Universal

INTENDED USE: Ideal for cementing crowns, bridges, inlays, onlays and posts made of materials such as: ceramic, metal, zirconia or resin.

INDICATIONS: In cases where maximum durability is required, aesthetic restorations, cases where adhesion is critical, and when easy handling is desired. When the light can't penetrate in the restoration **INTENDED**

USER: Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must make restorations with Medental Dual Cure Resin Cement.

PRECAUTIONS:

Some patients can present allergies to the components of the formula, in this case, suspend its use.

Keep away from Sunlight.

Keep away from rain.

Keep away from the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

No requires refrigeration.

Only for professional dental use.

Room temperature before use.

Not require cold chains to transport.

Do not use it after the date of expiry.

CONTRAINDICATIONS: The material must not be used when the application techniques described are not possible. Composite resins are subject to irritation in soft tissues and cells in direct contact when the monomers included in their formulation are in direct contact with these tissues, causing irritation and even pulp death if they are not 100% photocured. It is contraindicated to place composite resin in any presentation on materials that contain eugenol since these inhibit polymerization. Layers greater than 3 mm should not be placed to allow correct polymerization. The layers should not exceed 3mm to avoid fractures due to chewing compression. Humidity can affect adhesion. It should not be used in very thin or extra-thin spaces because it may change color. Self-etching materials can interfere with chemical self-curing, resulting in incomplete setting. In subgingival areas, it should be cleaned promptly to prevent gingival inflammation.

COMPOSITION: Crystalline and fumed silicon dioxide, Bis-GMA, Stabilizer and Photo initiators and chemical initiators.

INSTRUCTIONS FOR USE:

TECHNICAL INFORMATION: - Working time setting time may vary based on storage conditions, temperature, humidity, etc. - For better results it is recommended to use halogen or LED lamp with a wavelength between 400-500nm that activates the photoinitiators of the resin.

A. TOOTH PREPARATION: Follow usual procedures for tooth preparation. Exposed dentine must be

protected by a calcium hydroxide preparation. Glass ionomer or another suitable base can be used to cover the preparation. Remove temporary restorations isolate and clean the teeth of any traces of temporary cement.

B. ETCHING TOOTH SURFACE ENAMEL: Apply etchant over the surface enamel of tooth using a cotton pellet or brush. Allow to remain for 20-30 seconds. On primary tooth high on fluoride 90-120 seconds. If Glass Ionomer base is used, do not etch for more than 20-30 seconds. Flush area completely for 20 seconds. Do not allow patient to rinse the mouth or allow saliva to touch etched enamel surface. Air dry completely with an oil moisture free air source. The surface of tooth should have a dull white frosted appearance. If not, repeat the process allowing etchant to remain on tooth for an additional 20-30 seconds.

C. PRIMING TOOTH SURFACE: Apply a thin uniform coating of Enamel Bonding Resin immediately onto etched tooth surface. Thin the resin coating on the tooth with a gentle stream of moisture free air. Avoid contact with gingival or other soft tissue. When using a prime, apply a layer of porcelain primer on the etched porcelain veneers.

D.CEMENT MIXING AND APPLICATION: Dispense equal amount of Base and Catalyst on a mixing pad. Both base and catalyst should be immediately closed and kept in a cool, dark place. Mix the base and catalyst for 10 seconds to obtain a uniform mix. Coat the prepared teeth and restore completely with the mixture. Seat the restoration with a slightly rocking motion. Remove excess cement with a brush wetted with the bonding resin leaving a small strand of cement at the margin. Cure the cement with a wide diameter light for 40-60 seconds each of occlusal surface, proximal and gingival margins. Light transmitting wedges will help the gingival healing. Check the margins. Finish and polish the restoration carefully in appropriate manner. Adjust the occlusion, if necessary. NOTE: Translucent restorations (these which let through an adequate amount of light) need only be cemented by light-curing with the Base Cement.

CONSERVATION:

- The system is designed to be used at room temperature above 24°C or 75°F.
- Shelf life at room temperature is two years for unopened kits and one year for opened kit when syringe and vital caps are replaced immediately after dispensing.
- Unopened kits can be refrigerated (4°C-40°F) to extend shelf life. If stored in refrigerator, bring to room temperature prior to use.
- Do not store materials in proximity to eugenol containing products.
- Do not expose materials to elevated temperature or intense light.

WARNINGS and PRECAUTIONS: - In case of any change on the original function of the product, stop using it and contact the manufacturer. - The Dual-Cure Resin Cement is not irritable to the skin. - Read carefully and follow all instructions - Before the use of these products always use rubber dam gloves and protective glasses. - In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor. - In case of contact with skin wash the area with soap and water. - The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance. - Ask for patient to rinse his mouth immediately after treatment. - The product is designed to be used at room temperature. If refrigerated, the material may be difficult to extrude. - It is recommended to cover the polymerization time indicated by the manufacturer according to the layer thickness, because a non-polymerized resin can cause irritation or allergic reactions.

WARRANTY:

If the use of this product presents some anomaly such as:

- It does not fulfill the intended purpose,
- It doesn't return the functionality of the tooth,

-It doesn't have the properties it refers to,
-It has damage to packaging or manufacturing,
-The instructions for use are misunderstood,
-Any secondary reaction or anomaly:
Please send detailed report to: Conamco, S.A. de C.V.
email: atencionaclientes@conamcosacv.com.mx

Warranty: MEDENTAL and/or CONAMCO S.A. de C.V. warranty the product against any fault of manufacture of the material. MEDENTAL and/or CONAMCO S.A. de C.V. do not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective products.

The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at:

<https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS:

1. Use New and clean nitrile gloves.
2. The resin must be cured in a safe container for resin, either by exposing it to the sun or with a light cure lamp.
3. Once Cured, it can be discarded with solid residues.

NOTE: Liquid or partially cured resin must be disposed of in accordance with the regulations.

FEATURES/CLAIMS:

- Polymerization is assured.
- It is a biocompatible material if used correctly
- High mechanical resistance.
- Lower risk of failures due to lack of light.
- Excellent adhesion.
- Reduced postoperative sensitivity.
- Control of working time.
- Clinically versatile.

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- Veuillez envoyer un rapport détaillé à : Conamco, S.A. de C.V., email :