SELF ETCH DUAL CURE RESIN CEMENT

INSTRUCTIONS FOR USE

English
Español
Français

 CONAMCO, S.A. de C.V.
José Ma. Bustillos No. 28-Bis, 1st. and 2nd. Floor.
Col. Algarín, Alcaldía Cuauhtémoc.
P.C. 06880, Mexico City, Mexico.

 Ivonne Richaud.
Travessera de Gracia Street No. 315, 2nd.-1st.
P.C. 08025, Barcelona, Spain.

Rev."0" / 21 Jul 25

SELF ETCH DUAL CURE RESIN CEMENT

ENGLISH

PRODUCT DESCRIPTION

DESCRIPTION: Dual-cure resin is a dental material designed to polymerize by both light and chemical reaction, ensuring polymerization where light cannot fully penetrate.

CLINICAL BENEFITS:

- Working time
- Strong chemical adhesion
- Long-term duration.
- Easy to use (fewer clinical steps)

* Assured Polymerization

PERFORMANCE CHARACTERISTICS:

- High mechanical resistance
- Excellent adhesion.
- Long-term stability.

SHADES: Universal

INTENDED USE: Ideal for cementing crowns, bridges, inlays, onlays and posts made of materials such as: ceramic, metal, zirconia or resin.

INDICATIONS: In cases where maximum durability is required, aesthetic restorations, cases where adhesion is critical, and when easy handling is desired. When the light can't penetrate in the restoration **INTENDED USER:** Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must make restorations with Medental Self Eych Dual Cure Resin Cement.

PRECAUTIONS:

Some patients can present allergies to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away from the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

No requires refrigeration.

Only for professional dental use.

Room temperature before use.

Not require cold chains to transport.

Do not use it after the date of expiry.

CONTRAINDICATIONS: The material must not be used when the application techniques described are not possible. Composite resins are subject to irritation in soft tissues and cells in direct contact when the monomers included in their formulation are in direct contact with these tissues, causing irritation and even pulp death if they are not 100% photocured. It is contraindicated to place composite resin in any presentation on materials that contain eugenol since these inhibit polymerization. Layers greater than 3 mm should not be placed to allow correct polymerization. The layers should not exceed 3mm to avoid fractures due to chewing compression. Humidity can affect adhesion. It should not be used in very thin or extra-thin spaces because it may change color. Self-etching materials can interfere with chemical self-curing, resulting in incomplete setting. In subgingival areas, it should be cleaned promptly to prevent gingival inflammation.

COMPOSITION: Crystalline and fumed silicon dioxide, Bis-GMA, Stabilizer and Photo initiators and chemical initiators.

INSTRUCTIONS FOR USE:

A. APPLICATIONS:

1. Metal crowns, bridges, inlays and onlays (include porcelain fused to metal and composite to metal).
2. Resin crowns, bridges, inlays and onlays.
3. Metal (prefabricated or cast) and non-metal fiber endodontic posts.

4. Porcelain crowns, inlays and onlays (includes alumina and zirconia).
5. Veneer cementation.

B. TECHNICAL INFORMATION:

- **Working Time:** Minimum 1 min. 45 sec. (Includes mixing)
- **Setting time:** Maximum 4 min 30 sec.
- Working time setting time may vary based on storage conditions, temperature, humidity, etc.
- For better results is recommended to use halogen or LED lamp with a wavelength between 400-500nm that activates the photoinitiators of the resin.

C. INSTRUCTIONS FOR CROWNS, BRIDGES, INLAYS AND ONLAYS:

Preparation procedures:

NOTE: Resin Cement is self-adhesive cement which does not require etching or priming of the prepared surfaces.

1. Remove the temporary restoration and all temporary cement residue.
2. Try in the final restoration to insure proper fit.
3. Clean the preparation with pumice and water. Rinse thoroughly.
4. Remove all surface water using a strong stream of air on the preparation for 3 to 5 seconds. Do not desiccate. (Protect the preparation from contamination. When applicable, use a rubber dam is recommended).

Cement Mixing Procedures:

1. Remove the syringe cap or mixing tip from the dual syringe.

2. Dispense a small amount of material onto a mixing pad to eliminate any voids in each chamber of the dual syringe.
3. Attach a mix tip to the dual syringe by lining up the key and keyway. Then turn the brown locking mechanism clockwise.
4. Pressing the plunger will mix and dispense the Dual Cure Self-Etch.

Cementation Procedures:

NOTE: Since ceramics and porcelains vary greatly. It is important to consult with the ceramic/porcelain manufacturer or your laboratory for instructions on proper surface treatment.

1. Apply cement to the internal surface of the restoration. For inlays, it may be easier to apply the cement to the tooth preparation.
2. Seat the restoration and remove any excess cement.
HINT: To aid the removal of excess cement, initially light cure the margins for 2-3 seconds.
3. Light Cure for 20-30 seconds, or allow self-curing. (When using self-cured mode, light curing the margins even briefly will expedite the final set).

D. INSTRUCTIONS FOR POST (PRE-FABRICATED, METAL, CAST AND NON-METAL/FIBER POSTS).

Preparation Procedures:

1. Prepare the post space per post manufacturer's instructions.
Etching and bonding agents are not required.
2. Size and fit post to insure proper fit.
3. Clean post per manufacturer's instructions.

Cement Mixing Procedures:

1. Remove the syringe cap from the dual syringe.
2. Dispense a small amount of material onto a mixing pad to eliminate any voids in each chamber of the dual syringe.
3. Attach a mixing tip to the dual syringe by lining up the key and keyway. Then turn the brown locking mechanism clockwise.

Cementation Procedures:

1. Apply cement to the canal walls using an endodontic file, lentulo spiral or paper point and coat the post.
2. Seat the post and remove any excess cement.
HINT: To aid in the removal of excess cement, initially light cure the cement for 2-3 seconds.
3. Light cure for 20-30 seconds, or allow self-curing. (When using self-cured mode, light curing the margins even briefly will expedite the final set).
4. Proceed with core build-up procedure.

CONSERVATION:

- The system is designed to be used at room temperature above 24°C or 75°F.
- Shelf life at room temperature is two years for unopened kits and one year for opened kit when syringe and vital caps are replaced immediately after dispensing.
- Unopened kits can be refrigerated (4°C-40°F) to extend shelf life. If stored in refrigerator, bring to room temperature prior to use.
- Do not store materials in proximity to eugenol containing products.
- Do not expose materials to elevated temperature or intense light.

WARNINGS and PRECAUTIONS: - In case of any change on the original function of the product, stop using it and contact the manufacturer. - The Dual-Cure Resin Cement is not irritable to the skin. - Read carefully and follow all instructions - Before the use of these

products always use rubber dam gloves and protective glasses. - In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor. - In case of contact with skin wash the area with soap and water. - The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance. - Ask for patient to rinse his mouth immediately after treatment. - The product is designed to be used at room temperature. If refrigerated, the material may be difficult to extrude. - It is recommended to cover the polymerization time indicated by the manufacturer according to the layer thickness, because a non-polymerized resin can cause irritation or allergic reactions.

WARRANTY:

If the use of this product presents some anomaly such as:

- It does not fulfill the intended purpose,
- It doesn't return the functionality of the tooth,
- It doesn't have the properties it refers to,
- It has damage to packaging or manufacturing,
- The instructions for use are misunderstood,
- Any secondary reaction or anomaly:

Please send detailed report to: Conamco, S.A. de C.V.
email: atencionclientes@conamcosacv.com.mx

Warranty: MEDENTAL and/or CONAMCO S.A. de C.V. warranty the product against any fault of manufacture of the material. MEDENTAL and/or CONAMCO S.A. de C.V. do not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective products.

The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at:

<https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS:

1. Use New and clean nitrile gloves.

2. The resin must be cured in a safe container for resin, either by exposing it to the sun or with a light cure lamp.
3. Once Cured, it can be discarded with solid residues.

NOTE: Liquid or partially cured resin must be disposed of in accordance with the regulations.

FEATURES/CLAIMS:

- Polymerization is assured.
- It is a biocompatible material if used correctly
- High mechanical resistance.
- Lower risk of failures due to lack of light.
- Excellent adhesion.
- Reduced postoperative sensitivity.
- Control of working time.
- Clinically versatile.

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-Veuillez envoyer un rapport détaillé à : Conamco, S.A. de C.V., email :