RESTORATIVE SELF-CURE COMPOSITE

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RESTORATIVE SELF-CURE COMPOSITE

ENGLISH

PRODUCT DESCRIPTION: Polymerizes through a chemical reaction by mixing the Base and Catalyst in equal proportions.

CLINICAL BENEFITS:

- Can be used without a light-curing unit.
- Useful in pediatric dentistry and emergency procedures.

- Provisional restorations.

PERFORMANCE CHARACTERISTICS:

- Good or improved bond strength and adaptation after aging.
- Clinical alternative for some posterior restorations.
- Not uniformly equivalent to light-cure materials.
- Retention and short-term function.
- Medium shrinkage
- Working time: 2 min.

COLOR: Universal

INTENDED USE: Ideal for core build-ups and rapid Class I and Class II restorations. It is ideal in situations where a light-curing unit is not available.

INDICATIONS:

Fixation / cementation of restorations.
Fabrication of provisional restorations (crowns and bridges).
When the light polymerization is not possible.

INTENDED USER: Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

PATIENT POPULATION: All patients of any age, gender, weight, or ethnicity, with primary or permanent dentition, who, for different reasons, must be make restorations with Medental Restorative Self-Cure Composite.

PRECAUTIONS: Some patients may present allergy to the components of the formula, in this case, suspend it use.

Keep away from Sunlight.

Keep away from rain.

Keep away the reach of children.

Keep away from heat or sparks.

Temperature limit between 4°C and 24°C.

Fragile.

Non-Sterile.

Does not require refrigeration.

Only for professional dental use.

Room temperature before use. It does not require cold chain to transport.

Do not use after date of expiry.

CONTRAINDICATIONS: The material must not be used when the application techniques described are not possible. Composite resins are subject to irritation in soft tissues and cells in direct contact when the monomers included in their formulation are in direct contact with these tissues.

COMPOSITION: Crystalline and fumed silicon dioxide, Bis-GMA, Stabilizer and Chemical initiators.

INSTRUCTIONS FOR USE:

Composite filling material is a resin-bonded crystalline silicon dioxide composite. Starting with the original research of Dr. R.L. Bowen at the National Bureau of Standards, Washington, D.C., composite filling material has been developed to represent the finest composite dental restorative radiopaque material available. "Universal" filling material has superb optical properties.

A. CAVITY PREPARATION:

Prepare the cavity, remove the caries. If liner is required, composite filling material is compatible calcium hydroxide and Zinc Phosphate Cement. Etch the enamel surface in the apply layer of bonding prior the application of the composite.

B. ETCHING TOOTH SURFACE ENAMEL:

Apply etchant over the surface enamel of tooth using a cotton pellet or brush. Allow to remain for 20-30 seconds. On primary tooth high in fluoride 90-120 seconds. If Glass Ionomer base is used, do not etch for more than 20 30 seconds. Flush area completely for 20 seconds. Do not allow patient to rinse the mouth or permit saliva to touch etched enamel surface. Air dry completely with an oil/moisture free air source. The surface of tooth should have a dull, white frosted appearance. If not, repeat the process allowing etchant to remain on tooth for an additional 20-30 seconds.

C. MIXING INSTRUCTIONS:

1. Using the rounded end of the disposable plastic spatula, remove from the jar of Catalyst and amount of paste slightly less than volume of the cavity to be filled. Transfers this paste to the top sheet of the mixing pad. Replace cap on jar.
2. Using the flat squared end of the disposable plastic spatula, remove from the jar of Base an amount of paste approximately equal to the Catalyst Paste taken. Transfer this paste to the same sheet of the mixing pad alongside the Catalyst paste already there. Replace cap on jar.
3. Using the flat squared end of the disposable plastic spatula, mix the two pastes for about 20 seconds using 1 to 2 strokes per second. The objective is to make the mixture as homogeneous as possible during this time. The composite filling material is now ready to insert in the cavity.

D. MIXING PRECAUTIONS:

- Keep composite filling material jars covered when not in use.
- Remove only the amount of paste needed. Never return either composite filling material paste to its jar.
- Always use both ends of the spatula as described in "Mixing instructions".
- Be certain to avoid cross-contamination of the pastes. The presence of either paste in the other will result in hardening of the paste in the jars.
- Never use a disposable spatula a second time. A used spatula could be a source of contamination.
- Do not use metal spatula for mixing composite filling material.

E. PLACEMENT:

Insert the mixing composite filling material into the dried, prepared cavity using a non-metallic instrument designed for this purpose. Use a slight excess and apply matrix band. A mylar band is to be preferred composite filling material will begin to harden about 2 minutes from beginning of the mixing tier will allow about 1 minute 40 seconds for placement. The composite filling material will set during the next 2 minutes and should not be disturbed during this period. Another 2 minutes will complete the hardening time when the matrix band may be removed, and the surface of the filling finished and polished. See schematic "clock":

F. FINISHING: At the end of 6 minutes from beginning of mixing composite filling material. Remove feather flash with a sharp instrument. Finish if necessary, using

a diamond or carborundum stone or carbide bur various polishing materials are available which give good result.

CONSERVATION:

- The system is designed to be used at room temperature above 24°C or 75°F.
- Do not expose materials to elevated temperature or intense light.
- Close the jars when not in use.
- Once mixed, the product should not be stored.

WARNINGS and PRECAUTIONS: In case of any change on the original function of the product, stop using it and contact the manufacturer.

- The Restorative Self-Cure Composite is not irritable to the skin.
- Read carefully and follow all instructions
- Before the use of these products always use rubber dam gloves and protective glasses.
- In case of contact with the eye, wash immediately with plenty of water for at least 15 minutes and if the irritation persists consult a doctor.
- In case of contact with skin wash the area with soap and water.
- The product should not be ingested or aspirated. If accidentally swallowing induces vomiting immediately and get medical attention. In case of inhalation, clean your nose with water and seek medical assistance.
- Ask for patient to rinse his mouth immediately after treatment.
- The product is designed to be used at room temperature. If refrigerated, the material may be difficult to extrude.
- It is recommended to cover the polymerization time indicated by the manufacturer according to the layer thickness, because a non-polymerized adhesive can cause irritation or allergic reactions.

WARRANTY: If the use of this product on a patient or anyone else causes some anomaly, please send detailed report to: Conamco, SA de CV email:

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Warranty: MEDENTAL INTERNATIONAL and/or CONAMCO SA de CV warranty the product against any fault of manufacturing of the material. MEDENTAL INTERNATIONAL AND/OR CONAMCO SA de CV does not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product. The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at: <https://ec.europa.eu/tools/eudamed>

DISPOSAL CONSIDERATIONS: New, clean nitrile gloves should be worn.

1) *Removal of excess.*

- Fresh (before complete setting): Remove with a probe or spatula. recommends cleaning immediately.
- Setting: Using cement chisels, curettes or metal excavators with controlled movements, pass waxed dental floss in interproximal areas, pulling downwards to avoid fracturing restorations.

2) *Removal under restorations or crowns.*

- It is not possible to remove it without removing the restoration; in these cases, the crown must be sectioned with a carbide/diamond bur and the tooth cleaned manually or with ultrasound.

*In both cases, it is recommended to protect soft tissues, work with good lighting, and use a magnifying glass to identify cement residues.

3) *Dispose of with solid waste in accordance with the regulations of each country.*

FEATURES CLAIMS:

- Fracture resistance.
- Good adhesion to zirconia and metal.
- Emergency restorations.

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